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Instructions for Use TPP Tissue Culture Flasks with Re-closable Lid

TPP tissue culture flasks with re-closable lid are designed for the cell/tissue cultivation and growth in manual use. They meet the need for a large access from the above to the cell culture. The growth surface is opto-mechanically activated to optimize cell adhesion and proliferation. The angled flask neck reduces the risk of medium contact with the inside of the screw cap and allows excellent access for tissue culture scrapers and serological pipettes.

TPP cell culture flasks with re-closeable lid are available with filter screw caps only, when using a VENT screw cap a possible pressure build-up of 0.03 bar or higher could press open the re-closeable lid.

The TPP tissue culture flask with re-closeable lid is intended for single use.

Note:

Follow national regulations when handling biological material, use appropriate protective clothing

Take note of the regulations of aseptic working during the work process.

Instructions:

- *Check the packaging material and the product for intactness. Use only faultless material.*
- *Check the expiry date (EXP) on the label and packaging. Only use products with valid EXP.*
- *Open the packaging in a sterile environment and remove a complete system for use.*
- *Open the filter screw cap or the lid of the cell culture flask and fill the flask with medium and inoculum according to your standard laboratory routine.*
- *Close the filter screw cap completely. A continuous gas exchange takes place through the integrated hydrophobic membrane [0.22 µm]. If the PTFE filter membrane gets wet, the gas exchange is reduced for a short time.*
- *Cultivate and harvest cells according to your common laboratory routine.*

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
Opening Mechanism:

The lid can be opened and closed multiple times and remain gas- and liquid-tight.

- 1 Press here
- 2 Lift out lid

Important: To do this, it is necessary to place the flask on a firm surface. To open the lid, press the notch provided for this purpose with your forefinger. This allows you to control the pressure of your finger on the opening mechanism. Too much force can lead to breakage.

Advice and Recommendations:

- There are inscription fields on both sides of the flask for labelling.
- The graduation serves as a reference for the filling quantity and is not an absolute value.
- The stacking rim allows a stable stacking of several flasks on top of each other, also in combination with other TPP flasks of the same size.
- For a homogenous cultivation of the cells, levelling of the incubator and continuous supply of medium is essential.
- Avoid and/or eliminate vibration sources near the incubator to ensure uniform cell growth.

Note:

The barrier limits the growth area to 115 cm².

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Additional:

Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage www.tpp.ch.