Protocolo de Congelamento de Células Aderentes

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Cryopreservation of Adherent Cells

An essential prerequisite for the constant quality of cells is the establishment of a cell bank by means of cryopreservation. Continuous cultivation may lead to genetic changes, reduction of the proliferation rate, transformation, and altered expression patterns. With the help of cryopreservation, cell lines can be stored almost indefinitely without any reduction in quality.

The following protocol describes the basic steps for freezing adherent cells. This protocol is not suitable for every cell line. Please also refer to the protocols of the corresponding cell banks/sources.

Refer always to the national regulations for handling biological material, use the appropriate protective clothing and observe the rules of aseptic work in all steps.

General Information:

Before cells are frozen, they should be examined morphologically and checked for contamination. The cells should be in a late log phase. For adherent cells, cell harvesting at 70-80% confluence is recommended. 24 hours before freezing a media change should be performed.

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Freezing Media:

The following media are suitable as freezing media:

Media with Serum:

- 90% culture media +10% DMSO*
- 45% sterile filtered conditioned media / 45% fresh media + 10% DMSO*.

The antifreeze in the culture media, and the high protein content of the serum, protect the cells from damage during the freezing and thawing process. With a serum-free media the protective effect of the serum is lost and can be compensated by additional antifreeze (e.g.: HES, Pluronic F68™, and methyl cellulose).

Media Serum-free: (xeno-free, chemically defined)

- 89% chemically defined, serum-free media, 10% DMSO*, 1% Pluronic F68™ [5]
- Commercially available serum-free freezing media

**DMSO is not suitable for some cells HL-60 [1], HBE4-E5/E7 [1, 2], glycerol is be a suitable alternative*

Preparation:

Pre-cool the serological pipettes, the freezing media, the freezing container and the cryotubes in the refrigerator.

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Procedure:

- *Remove all the culture media from the flask and wash the cells twice with warm PBS w/o Ca and Mg (0.1 mL/cm²).*
- *Dissolve the cells with a suitable dissociation reagent (e.g. Trypsin, Trypsin-EDTA). Follow the protocol established in your laboratory.*
- *After addition of fresh media containing a dissociation inhibiting component, transfer the cells into a centrifuge tube.*
- *Resuspend the cells with fresh medium and determine the total cell count and live cell count (viability).*
- *Calculate the required volume of freezing medium to achieve a final live cell count of 1×10^6 - 1×10^7 cells/mL*
- *Centrifuge the cells at 200 x g for 5 min.*
- *Remove the supernatant and resuspend the cells with the calculated volume of freezing medium.*
- *Aliquot the cell suspension into the appropriately labelled cryotubes.*
- *The cells should be frozen at a controlled freezing rate of -1 °C/minute. You can use programmable freezers or freezing container.*
- *Store the freezing container overnight in a -80 °C freezer.*
- *Transfer the cryotubes the next day to a -152 °C freezer or store them in the gas phase of an LN₂ cryogenic storage unit. Since no reactions take place in the cells from -130 °C [7], storage in the gas phase is possible and recommended for safety reasons [8]*

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Thawing of the Cells:

- *Fill 14 mL warm culture media without antibiotics in a T-75 cm² flask and 10 mL in a centrifuge tube.*
- *Remove one cryotube from the LN2 cryogenic storage unit and warm it up in a warm water bath equilibrated to the culture temperature of the cell line.*
- *As soon as ice crystals are no longer visible, transfer the content of the cryotube into the centrifuge tube using a sterile pipette.*
- *Centrifuge the cells: 200 x g for 5 minutes.*
- *Remove the supernatant and resuspend the cells with 1 mL fresh media.*
- *Avoid foaming during resuspension. This can lead to poor and uneven cell attachment. Foaming can be reduced by an even pipetting rhythm which is not performed quickly and hastily.*
- *Transfer the cells into the T-75 cm² flask and cultivate for 24 hours.*
- *Transfer the cells into the T-75 cm² flask and cultivate for 24 hours.*
- *After 24 hours you can check the cells for vitality, identity, and sterility and sub-cultivate them if desired.*
- *The thawing controls should be performed regularly.*

1. ATCC® Animal Cell Culture Guide. Tips and Techniques for Continuous Cell Lines: 24-27, 2014
2. ECACC Protocol Fundamental Techniques in Cell Culture. Laboratory Handbook 2nd Edition: 51-52, 2010
3. Freshney R., Culture of Animal Cells: A Manual of Basic Technique and Specialized Applications. Wiley-Blackwell: 307-326, 2016
4. Gstraunthaler G., Lindl T., Zell- und Gewebekultur: Allgemeine Grundlagen und spezielle Anwendungen. Springer Spektrum Verlag: 137-140, 2014
5. Hernández, Y.G., Fischer, R. Serum-free Culturing of Mammalian Cells – Adaptation to and Cryopreservation in Fully Defined Media. ALTEX 24: 110-116, 2007
6. Morris C.B., Cryopreservation of Animal and Human Cell Lines. In: Day J.G., Stacey G.N. (eds) Cryopreservation and Freeze-Drying Protocols. Methods in Molecular Biology™, vol 368. Humana Press 2007
7. Stacey G., Masters J., Cryopreservation and banking of mammalian cell lines. *Nat Protoc* **3**, 1981–1989 2008
8. Tedder R., Zuckerman M., Hepatitis B transmission from contaminated cryopreservation tank. *The Lancet* 346(8968), 137–140. 1995