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---

### TubeSpin® Bioreactor 50 ml, TPP product number 87050

1. **Take a new sterile TubeSpin® Bioreactor 50 ml**  
*Each experiment is done in a new sterile bioreactor. The geometry of the individual bioreactors is absolutely identical.*
2. **Dilute cells into pre-warmed medium and add 5 ml per bioreactor**  
*5 ml is suggested for routine applications, but other volumes work as well.*
3. **Agitate the bioreactor at 150 rpm in a shaker that provides 5% CO<sub>2</sub>**  
*In the presence of CO<sub>2</sub> and humidity (at least 80 %) the 5 filter openings (A, B, C, D, E) stay "open". If no CO<sub>2</sub> control is available, the tubes can be easily "closed" by covering the openings with an adhesive tape. If the pH starts to drop, the tape can be removed from position "A". Closing the openings "B" to "D" also helps to limit evaporation (e.g. no humidity control, small working volumes). The opening „A“ is designed to provide sufficient gas exchange for such cases.*
4. **Take samples to count cells and analyze the culture**  
*The tubes can be sampled the same way as any other culture vessel. Alternatively, several identical tubes are prepared and the whole tube is sacrificed as a sample. This strategy has the advantage that the sampling procedure has no impact on the remaining tubes.*
5. **Harvest of cells or supernatant by centrifugation**  
*Centrifugation of the bioreactor is convenient and provides cells and/or supernatants without any risk of contaminations.*