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Instruction for Use TPP PCV tube (Packed Cell Volume)

TPP „easy read“ Measuring Device

The PCV tube, together with the "easy read" measuring device, allows easy and fast measurement of biomass in cell culture, microbiology and suspension cultures of fungi and yeasts. You need a centrifuge with swing-out rotor.

The PCV tube is intended for single use, the "easy read" measuring device is reusable.

Note:

Follow national regulations when handling biological material, use appropriate protective clothing

Take note of the regulations of aseptic working during the work process.


Instructions PCV tube

- *Check the packaging material and the product for intactness. Use only faultless material.*
- *Use only PCV tubes with a valid expiration date (EXP).*
- *Transfer a sample of 100 - 1000 µl of a mixed cell suspension into the PCV cell measuring tube.*
- *Close the tube with the cap #87008 if necessary.*
- *Centrifuge the tube for 1 minute at 2'500 x g in a swing-out rotor for 1.5 | 2 mL vessels.*
- *Read the volume of the cell pellet as described in the TPP "easy read" measuring instrument instructions.*
- *Discard the PCV tube*

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Instruction "easy read" Measuring Device

- After centrifugation, insert the PCV tube ① into the opening provided in the "easy read" after centrifugation.
- Move the magnifying glass ③ until the height of the cell pellet in the capillary is equal to the cross.
- Read the cell pellet volume in μl with the "easy read" measuring device.

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Cleaning and disinfection of the "easy read"

Please note the following when choosing the cleaning and disinfecting agent:

- Choose procedures with neutral, aluminum-compatible cleaning and disinfecting agents and fully demineralized water.
- Avoid abrasive cleaning agents or scratching metal brushes, agents containing halogen or chloride: risk of corrosion.
- Thermal cleaning (autoclaving) is not recommended.
- Remove all cleaning agent residues by rinsing thoroughly.
- After cleaning/disinfection, dry the instrument completely and store it in a dry place.

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x = Densidade celular número de célula/mL

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Comments and recommendations:

- The PCV method does not provide information about the cell viability of the culture.
- A new calibration line must be determined for each cell type.
- Cell pellet cannot be removed from the capillary.

Calculation PCV (%) value:

- The measured volume is expressed in % of the volume of the transferred suspension sample. It corresponds to the total biomass of the culture.

$$PCV(\%) = \frac{\text{volume cellpellett}}{\text{volume probe}} \times 100$$

Calculation of cell density:

To determine the cell density from the biomass (PCV %) a calibration line must be established.

1. Take a 2 ml sample from a well mixed cell suspension.
2. Carry out a double determination for the PCV value (%) and for the cell density using the counting chamber (cell number/ml).
3. The frequency of sampling and measurements performed depends on the growth rate (doubling time) of the cells.
4. Enter the determined mean values in a table and create a diagram (see Fig. 2).
5. Determine the linear equation based on the measured values (Fig. 2)

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
7. The cell density of a sample cell number/ml can be calculated from the determined straight line equation:

$$\text{Cell density} = \frac{(\text{PCV} - b)}{m}$$

Fig.2: A typical calibration line for a cell line (NS0) which correlates the cell count per millilitre (determined by microscopic counting) with the biomass value of the sample expressed in percent.

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Materials

Measuring device aluminum
Magnifying lens glass

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Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage www.tpp.ch

Literatura:

1. Stettler, M., Jaccard, N., Hacker, D., De Jesus, M. Wurm, F. M., Jordan, M. (2006): New disposable tubes for rapid and precise biomass assessment for suspension cultures of mammalian cells. *Biotechnol. Bioeng.*95 (6): 1228–1233.
2. Wurm, F., Tanner, R., Jordan, M. (2008) Eine neue Methode zur Bestimmung von Biomasse in der Zellkultur. *BioSpektrum* 05.08 508-509