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## Instructions for Use TPP Tissue Culture Flasks

*TPP tissue culture flasks are intended for tissue culture in manual use. For optimal cell adhesion and growth, only the bottom of the tissue culture flask is opto-mechanically activated. The angled neck reduces the risk of medium contacting the inside of the screw cap by spillage and allows excellent access for cell scrapers and serological pipettes.*

*TPP tissue culture flasks are available with filter or VENT screw caps.*

*The tissue culture flask is for single use only.*

### Note:

*Follow the national regulations for handling biological materials and wear the appropriate protective clothing.*

*During the work process, be sure to observe the rules of aseptic technique.*

### Instructions:

- *Inspect packaging and product for integrity. Use only materials in perfect condition.*
- *Check the expiry date (EXP) on the label and packaging. Use only products with a valid EXP.*
- *Open the package in a sterile environment and remove a complete system for use*
- *Open the flask and fill it with the medium and inoculum according to your laboratory routine. Please refer to the optimal filling volume (see technical data).*

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**Filter Screw Cap:**

- Close the filter screw cap with ventilation holes. Continuous gas exchange occurs through the integrated 0.22 µm hydrophobic membrane (Figure 1). If the membrane becomes wet, gas exchange is temporarily reduced.


Figure 1: Filter-screw-cap closed tightly

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**VENT Screw Cap:**

- For a gas exchange, turn the VENT cap into the ventilation position at 12 o'clock. You will hear/notice a click.

- Visual check:  
Rectangle at 12 o'clock (Position 1)

During incubation of cells, the rectangle must be in position 1 / 12 o'clock.


- To stop gas exchange, rotate the cap 1/4 turn clockwise beyond the vent position. The flask is now tightly sealed and there is no gas exchange.

- Visual check:  
Rectangle at 3 o'clock. (Position 2)

- For a transfer of the flask within in the laboratory (laminar flow ↔ incubator), TPP recommends closing the VENT screw cap tightly.

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
#### **Advice and Recommendations:**

- *For labelling there are bilateral inscription fields*
- *The volume scale is an indication of the filling volume, not an absolute value.*
- *The stacking rim allows a stable stacking of several flasks on top of each other, also in combination with other TPP flasks of the same size.*

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**Additional:**

*Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from [www.tpp.ch](http://www.tpp.ch).*