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TPP PCV-Röhrchen

(Packed Cell Volume)

TPP „easy read“ Messinstrument

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TPP PCV (Packed Cell Volume) tube

TPP “easy read” Measuring Device

The PCV tube, together with the “easy read” measuring device, allows easy and fast measurement of biomass in cell culture, microbiology and suspension cultures of fungi and yeasts. You need a centrifuge with swing-out rotor.

The PCV tube is for single use only. The “easy read” measuring device is reusable.

Note:

Follow the national regulations for handling biological materials and wear the appropriate protective clothing.

During the work process, be sure to observe the rules of aseptic technique.

Instructions TPP PCV tube

- *Inspect packaging and product for integrity. Use only materials in perfect condition.*
- *Check the expiry date (EXP) on the label and packaging. Use only products with a valid EXP.*
- *Transfer a sample of 100 - 1000 µL of a mixed cell suspension into the PCV cell measuring tube.*
- *Close the tube with the cap # 87008 if necessary.*
- *Centrifuge the tube for 1 minute at 2'500 x g in a swing-out rotor for 1.5 | 2 mL vessels.*
- *Read the volume of the cell pellet as described in the TPP “easy read” measuring instrument instructions.*
- *Discard the PCV tube.*

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
Instruction "easy read" Measuring Device

- After centrifugation, insert the PCV tube ① into the opening provided in the "easy read" after centrifugation.
- Move the magnifying glass ③ until the height of the cell pellet in the capillary is equal to the cross.
- Read the cell pellet volume in µL with the "easy read" measuring device.

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Calculation PCV (%) value:

- The measured volume is expressed in % of the volume of the transferred suspension sample. It corresponds to the total biomass of the culture.

$$\text{PCV}(\%) = \frac{\text{volume cellpellett}}{\text{volume probe}} \times 100$$

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Cleaning and disinfection:

“easy read“-body:

- Wipe with a soft cloth dampened with a small amount of neutral detergent to remove stains.
Then wipe the “easy read” clean with a cloth dampened with lukewarm water. Be careful not to touch the magnifier when cleaning.
- Avoid abrasive cleaning agents, scratching metal brushes, agents containing halogens or chlorides: risk of corrosion.
- For disinfection, TPP recommends 70% ethanol.

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“easy read”- magnifying glass:

- *Use compressed air to remove dirt and dust adhering to the surface.*
- *Then take a piece of lens paper and fold it around your finger to form a tip. Be sure to use lens paper for this step, not lab tissues or paper towels. Traditional tissues contain loose, coarse fibers that can scratch or leave marks on the lens.*
- *Put a small amount of 70% ethanol on the tip of the lens paper.*
- *When cleaning the lens surface, wipe in a circular motion from the inside to the outside.*
- *Allow the “easy read” to dry completely after cleaning/disinfection. TPP recommends dust-free and dry storage of the device.*

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Comments and recommendations:

- *Sample dilution is not required even at high cell densities. A reduction of the selected sample volume is sufficient.*
- *The PCV method does not provide direct information on the cell viability of the culture.*
- *Cell pellet cannot be removed from capillary*

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Technical Data:

Materials

Measuring device aluminum
 Magnifying lens glass

Abmessungen	Measurements	87010
Volumengraduierung μl	Volume graduation μL	0.0 – 5.0
Länge x Breite x Höhe mm	Length x width x height mm	320 x 15 x 13.5
Lupe: Länge x Breite x Höhe mm	Magnifying lens: Length x width x height mm	30 x 33 x 38

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Additional:

Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from www.tpp.ch.

Literature:

1. Stettler, M., Jaccard, N., Hacker, D., De Jesus, M. Wurm, F. M., Jordan, M. (2006): New disposable tubes for rapid and precise biomass assessment for suspension cultures of mammalian cells. *Biotechnol. Bioeng.* 95 (6): 1228–1233.
2. Wurm, F., Tanner, R., Jordan, M. (2008) Eine neue Methode zur Bestimmung von Biomasse in der Zellkultur. *BioSpektrum* 05.08 508-509