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Instructions for Use TPP Syringe Filter

TPP sterile syringe filters 22 and 45 are designed for filtration of aqueous solutions such as cell culture media and biological fluids from microorganisms, particles and undissolved precipitates larger than 0.22 µm or 0.45 µm - depending on the filter. They can be used bidirectional: Forward (syringe to container) or reverse (container to syringe), but only in one direction per filter unit and application.

The syringe filter is for single use only.

Note:

Follow the national regulations for handling biological materials and wear the appropriate protective clothing.

During the work process, be sure to observe the rules of aseptic technique.

Instructions:

- *Inspect packaging and product for integrity. Use only materials in perfect condition.*
- *Check the expiry date (EXP) on the label and packaging. Use only products with a valid EXP.*
- *Open the package in a sterile environment and remove the filter for use.*

Forward Filtration:

- *Take a syringe and remove the syringe plunger.*
- *Under sterile conditions, remove the backing from syringe blister pack.*
- *Place the syringe on the Luer-Lock or cone fitting of the syringe filter, grasping it as you place it on the blister.*

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- *Fill the syringe with the solution to be filtered.*
- *Carefully reinsert the syringe plunger into the syringe body and filter the liquid by gently pressing the plunger into an appropriate sterile container.*

Reverse filtration:

- *Remove a syringe from the blister pack.*
- *First, draw some air into the syringe to prevent a possible backflow of the liquid, then draw the liquid to be filtered from the container into the syringe.*
- *Remove any adhering liquid residue from the syringe tip.*
- *Under sterile conditions, remove the backing from syringe blister pack.*
- *Place the syringe on the Luer-Lock or cone fitting of the syringe filter, grasping it as you place it on the blister*
- *Filter the liquid by evenly pressing the syringe plunger into a suitable sterile container.*

Remarks and Recommendations:

- *The diameter of the TPP syringe filter fits perfectly on a standard 50 mL centrifuge tube. (Figure 1)*

Figure 1

- *Do not use the TPP filter at temperatures above 45 °C.*

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Technical Data:

Materials

Membrane Polyethersulfone (PES),
 type "TPP fast flow"
 Housing Polyethylene terephthalate (PET)

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Chemical Compatibility:

Compatible with most aqueous solutions, see
 documentation at www.tpp.ch.

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Additional:

Instructions for use, chemical resistance lists,
 and quality certificates for each product are
 available for download from www.tpp.ch.