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DiHEMDA	di(2-hydroxyethyl)methyldodecylammonium
DMSO	Dimethylsulfoxid
HPLC	High Performance Liquid Chromatography
HRMS	High Resolution Mass Spectrometry
LC	Liquid Chromatography
MeOH	Methanol
n.a.	Not Applicable, analyte not found
ng	Nanogram
NTS	Non-Target-Screening
NVOC	Non-Volatile Compound
USP	United States Pharmacopeia
UV	Ultra Violet

Literatur

- [1] McDonald, G Reid et al. "Bioactive contaminants leach from disposable laboratory plasticware." *Science (New York, N.Y.)* vol. 322,5903 (2008): 917. doi:10.1126/science.1162395
- [2] McDonald, Reid et al. "Bioactive Leachates from Lab Plastics Use of Plastic Disposables May Compromise Bioassay Results." *GIT laboratory journal Europe* 13 (2009): 24-26.
- [3] Olivieri, Aldo et al. "On the disruption of biochemical and biological assays by chemicals leaching from disposable laboratory plasticware." *Canadian journal of physiology and pharmacology* vol. 90,6 (2012): 697-703
- [4] Watson, John et al. "Extraction, identification, and functional characterization of a bioactive substance from automated compound-handling plastic tips." *Journal of biomolecular screening* vol. 14,5 (2009): 566-72.
- [5] Lee, Tet Woo et al. "Chemicals eluting from disposable plastic syringes and syringe filters alter neurite growth, axogenesis and the microtubule cytoskeleton in cultured hippocampal neurons." *Journal of neurochemistry* vol. 133,1 (2015): 53-65.
- [6] Schauer KL, Broccardo CJ, Webb KM, Covey PA, Prenni JE. Mass spectrometry contamination from Tinuvin 770, a common additive in laboratory plastics. *J Biomol Tech.* 2013 Jul;24(2):57-61
- [7] Lewis, L Kevin et al. "Interference with spectrophotometric analysis of nucleic acids and proteins by leaching of chemicals from plastic tubes." *BioTechniques* vol. 48,4 (2010): 297-302.
- [8] Budde, Dana et al. "Interaction of leachable model compounds and their impact on Chinese hamster ovary cell cultivation." *Biotechnology progress* vol. 37,4 (2021)
- [9] Eibl R, Steiger N, Fritz C, et al. Recommendation for leachables studies: standardized cell culture test for the early identification of critical films for CHO cell lines in chemically defined culture media. *Dechema.* 2014;1-22.

Disclaimer

TPP products are for research use only and not for clinical, diagnostic or therapeutic use. All products are intended for use by trained personnel that are familiar with safe laboratory practices.

TPP assumes no responsibility for damage or defects resulting from improper or unauthorized use. It is the responsibility of the user to store, handle, and use the products in accordance with the instructions provided.

TPP does not warrant the completeness or accuracy of this TechDoc. TPP's recommendations are intended as general guidelines and may not cover all possible scenarios. TPP shall not be liable for any indirect, incidental, consequential, or special damages arising out of the use or inability to use the information in this TechDoc.

Swiss law governs these terms of use and any resulting legal matters.