




Instructions for Use in Pictures TPP Syringe Filter



Do not re-sterilize or reuse the TPP syringe filter.
 For single use only.

Measurements	99722	99745
Inlet to outlet mm	27	27
Diameter mm	33	33
Filter size cm ²	4.52	4.52
Pore size µm	0.22	0.45
Max. using temp °C	45	45
Flow rate water 25 °C / 2 bar	≥ 150 mL/min	≥ 180 mL/min
Burst pressure at 25 °C	10 bar	10 bar
Filtration volume mL	10 – 200	10 – 200
Dead volume following further flushing with air mL	0.1	0.1
Connection In-let	Luer-lock	Luer-lock
Connection out-let	Luer-conus	Luer-conus