





















Recycling Codes

Note: Used products contaminated with hazardous waste must not be placed in the recycling bin, as they pose health risk to workers and the facility.

The ASTM International Resin Identification Coding System (RIC) is used to classify and sort products or packaging by material type, in accordance with EU Directive 97/129/EC.

Code	Product Raw Material	Number
	Polyethylene terephthalate (PET (filter membrane))	01
	High-density polyethylene (PE-HD)	02
	Low-density polyethylene (PE-LD (bags and films))	04
	Polypropylene (PP)	05
	Polystyrene (PS)	06
	All other plastics e.g., PTFE (membrane filter screw cap)	07
	Lithium battery	12
	Cardboard	20
	Paper	22
	Aluminum (easy read, adapter tray)	41
	Wood	50



Product Packaging per Product-Group		
Code	Product Raw Material	Number
	Pouch centrifuge tubes, racks	04
	Stretch film, 23 my (pallet wrapping)	04
	Blister plastic / plastic (pipette)	04
	Mixed material, not separable into individual materials Blister syringe filter Blister Turbo-Fix filter Blister Clipmax / Slidemax Blister paper / plastic (scraper, lifter, pipette) Peel-open-bags: 3-B Triple bag products Cryo tube Filtermax PCV-tube Tissue culture dish Tissue culture flask Tissue culture test plate Tissue culture tube	7
	TPP cardboard boxes Turbo-Fix	20
	Packaging material paper	22
	Wooden pallet	50

Disclaimer

TPP products are intended for Research Use Only (RUO) and are not approved for clinical, diagnostic, or in vitro fertilization (IVF) applications. The full Terms & Conditions, including limitations of warranty and liability, intended use, and reseller obligations, are available at:

https://www.tpp.ch/page/qualitaets_sicherung/index.php

Distributors who purchase and distribute TPP products acknowledge and agree to these Terms & Conditions and the associated disclaimer.